



Discovery through partnership | Excellence through quality

**StressXpress®**

## **20S Proteasome Activity Kit**

Catalog# SKT-133 (96-Well Kit)

Fluorometric detection of purified 20S proteasome activity

# TABLE OF CONTENTS

<b>GENERAL INFORMATION</b>	3	Materials Supplied
	4	Precautions
	4	If You Have Problems
	4	Storage and Stability
	5	Materials Needed but Not Supplied
<b>INTRODUCTION</b>	6	Background
	7	About this Assay
<b>PRE-ASSAY PREPARATION</b>	8	Overview and Assay Preparation
<b>ASSAY PROTOCOL</b>	12	Assay Protocol
<b>ANALYSIS</b>	14	Data Analysis
	16	AMC Standard Curve
<b>RESOURCES</b>	18	Plate Configuration Table
	19	References
	20	Warranty and Limitation of Remedy
	21	Notes

## GENERAL INFORMATION

### Materials Supplied

Catalog Number	Reagent	Quantity
SKC-133A	10x StressXpress® proteasome assay buffer	1 mL (reconstituted)
SKC-133B	20S proteasome	10µg
SKC-133C	MG132 (control inhibitor)	100µg
SKC-133D	StressXpress® 96 well assay plate	1 plate
SKC-133E	StressXpress® AMC standard	1.75µg
SKC-133F	Suc-LLVY-AMC (chymotrypsin-like substrate)	500µg

If any of the items listed above are damaged or missing, please contact our Customer Service department at (250) 294-9065. We cannot accept any returns without prior authorization.



**WARNING: Not for human or animal disease diagnosis or therapeutic drug use.**

## Precautions

**Please read these instructions carefully before beginning this assay.**

Kit contains sufficient materials for one 96-well plate assay or up to 100 x 50 $\mu$ L micro-cuvette based assays, controls and optional calibration standards. The reagents in this kit have been tested and formulated to work exclusively with StressMarq Biosciences Inc.'s StressXpress<sup>®</sup> Kits. This kit may not perform as described if any reagent or procedure is replaced or modified.

**For research use only. Not for human or diagnostic use.**

## If You Have Problems

### **Technical Service Contact Information**

<b>Phone:</b>	250-294-9065
<b>Fax:</b>	250-294-9025
<b>E-Mail:</b>	techsupport@stressmarq.com
<b>Hours:</b>	M-F 9:00 AM to 5:00 PM PST

In order for our staff to assist you quickly and efficiently, please be ready to supply the lot number of the kit (found on the outside of the box).

## Storage and Stability

The StressXpress<sup>®</sup> 20S proteasome assay kit can be stored as supplied at 4°C for up to 3 months, until ready to use, to ensure stability and activity. For storage of kit component solutions see 'Assay preparation'.

## Materials Needed But Not Supplied

- Fluorescence microplate reader with filters for 360nm excitation and 460nm emission
- Microfuge tubes (1.5mL)
- Sterile water (dH<sub>2</sub>O)
- Dimethyl sulfoxide (DMSO)

## Background

Proteasomes are non-lysosomal proteolytic complexes localised primarily in the cytoplasm and in the nucleus of eukaryotic cells. They are responsible for the ubiquitin-mediated degradation of short half-life proteins and peptides that are involved in essential cellular processes including cell-cycle regulation, apoptosis and transcriptional regulation, innate immunity and antigen processing, and in the removal of redundant or damaged proteins. As such protein degradation by the ubiquitin-proteasome pathway has a major regulatory function for proliferation activity and survival of both normal and malignant cells, and its dysfunction has been implicated in a wide range of other disease processes including neurodegenerative, cardiovascular and metabolic disorders.

The 26S proteasome structure is composed of a 20S proteasome catalytic core complex and one or two 19S (PA700) regulatory subcomplexes. The 20S core comprises two copies of 14 subunits (7 alpha subunits and 7 beta subunits) arranged in a  $\alpha_7\beta_7\beta_7\alpha_7$  cylindrical array. Proteolytic activities are determined by the  $\beta_1$  (caspase-like),  $\beta_2$  (trypsin-like) and  $\beta_5$  (chymotrypsin-like) subunits, access to which is guarded by the  $\alpha$ -subunits. The 19S regulatory unit consists of six ATPase and at least ten non-ATPase subunits that are required for ubiquitinated protein binding, deubiquitination, substrate unfolding and translocation to the 20S catalytic core.

Varying catalytic subunit composition ( $\beta_1$ ,  $\beta_{1i}$ ;  $\beta_2$ ,  $\beta_{2i}$ ;  $\beta_5$ ,  $\beta_{5i}$ ) results in a variety of possible subtypes from full constitutive proteasome ( $\beta_1$ ,  $\beta_2$ ,  $\beta_5$ ) through mixed populations to full inducible / immunoproteasome ( $\beta_{1i}$ ,  $\beta_{2i}$ ,  $\beta_{5i}$ ). Alternative regulatory complexes such as the PA200 and 11S proteasome activators confer different substrate specificities and activity compared to the 19S regulator.

## About This Assay

Flexible, robust, easy-to-use kit for measurement of 20S proteasome activity

- Measure 20S proteasome activity in continuous kinetic and end-point assays
- Utilise high sensitivity fluorogenic substrate for analysis of proteasome chymotrypsin-like ( $\beta 5$ ) activity
- Test 20S proteasome activators and inhibitors in micro-cuvette or 96-well plate HTS compatible format
- Quantify results using calibration standards supplied

The StressXpress® 20S proteasome assay kit facilitates the rapid, robust measurement of proteasome activity *in vitro*. The kit utilises high purity, fluorogenic substrate Suc-LLVY-AMC together with suitable calibration standards and controls for the accurate and sensitive assessment of proteasome chymotrypsin-like ( $\beta 5$ ) activity. Continuous kinetic or end-point assays can be performed in micro-cuvettes or in 96-well plate format for multi-sample analysis. Contains sufficient materials for one full 96-well plate assay or up to 100x 50 $\mu$ L micro-cuvette based assays to be run.

## Use this kit to

1. Screen potential proteasome inhibitors and activators
2. Assess proteasome chymotrypsin-like ( $\beta 5$ ) activity allowing relevant kinetic data (e.g.  $K_M$ ,  $k_{cat}$ ,  $k_{cat}/K_M$ ) or specific activity to be calculated
3. Evaluate performance of proteasome activators and inhibitors enabling relevant kinetic data to be determined (e.g.  $IC_{50}$  /  $EC_{50}$ , % inhibition / activation,  $K_i$ )
4. Investigate effect of activator complexes and binding proteins on proteasome activity
5. Optimise assays for specific user provided proteasome preparations to facilitate their use in high throughput screening (HTS)

## PRE-ASSAY PREPARATION

### Overview

The StressXpress® 20S proteasome assay protocol described can be utilised or adapted to:

- Assess performance of 20S proteasome, activators and inhibitors of interest enabling relevant kinetic data to be generated
- Screen potential proteasome inhibitors and activators
- Investigate effect of activator complexes and binding proteins on proteasome activity
- Run assays in micro-cuvette based format
- Optimise assays for specific user provided proteasome preparations to facilitate their use in high throughput screening (HTS)

### Assay Preparation

#### Proteasome assay buffer

Dissolve lyophilised 10× StressXpress® proteasome assay buffer supplied in 1mL dH<sub>2</sub>O and mix well by vortexing. Dilute resulting 10× StressXpress® proteasome assay buffer solution 10 fold with dH<sub>2</sub>O as required to give 1× proteasome assay buffer for use in assays and preparation of diluted solutions.

**Storage:** Proteasome assay buffer solutions can be stored short term at 4°C and longer term at -20°C or -80°C.



## Suc-LLVY-AMC (chymotrypsin-like substrate)

Recommend 100 $\mu$ M substrate assay concentration

- Suc-LLVY-AMC substrate (500 $\mu$ g) supplied lyophilised
- Prepare a 10mM stock solution by dissolving solid in 65.4 $\mu$ L DMSO
- Prepare a 1mM working solution by 10 $\times$  dilution in 1 $\times$  proteasome assay buffer – add 10 $\mu$ L 10mM Suc-LLVY-AMC stock solution to 90 $\mu$ L 1 $\times$  proteasome assay buffer to give 100 $\mu$ L working solution

**Storage:** Stock and working solutions can be stored at -80°C for up to three months protected from light. Avoid multiple freeze-thaws

## 20S proteasome

Recommend ~2.5nM 20S proteasome assay concentration (0.1 $\mu$ g / 50 $\mu$ L reaction)

- 20S proteasome supplied at 1mg/mL in 1 $\times$  proteasome assay buffer
- Prepare the required amount of 20 $\mu$ g/mL (~25nM) 20S proteasome by 50 $\times$  dilution with 1 $\times$  proteasome assay buffer – for example, add 10 $\mu$ L 1mg/mL 20S proteasome to 490 $\mu$ L 1 $\times$  proteasome assay buffer to give 0.5mL 20 $\mu$ g/mL 20S proteasome solution (5 $\mu$ L required per reaction)

**Storage:** DILUTED 20S PROTEASOME SOLUTION MUST BE PREPARED AND USED FRESH – DO NOT STORE. 20S proteasome may be stored as supplied for up to six months at 4°C. Recommend storage at 4°C if material is to be used on a regular basis. For longer term storage prepare single use aliquots, add high purity glycerol to 50% and store at -80°C. Avoid freeze-thaw cycles. Use frozen 20S proteasome aliquots immediately upon thawing to ensure activity retained.

**NOTE:** USER PREPARED 20S PROTEASOME CAN BE ANALYSED USING THE STRESSXPRESS® 20S PROTEASOME ASSAY KIT BY DILUTING IN 1 $\times$  PROTEASOME ASSAY BUFFER AND USING 2.5nM ASSAY CONCENTRATION AS A STARTING POINT

## MG132 (control inhibitor)

Recommend 5 $\mu$ M MG132 inhibitor assay concentration

- MG132 (100 $\mu$ g) supplied lyophilised
- Prepare a 500 $\mu$ M stock solution by dissolving solid in 420 $\mu$ L DMSO
- Prepare a 50 $\mu$ M working solution by 10 $\times$  dilution in 1 $\times$  proteasome assay buffer – add 10 $\mu$ L 500 $\mu$ M MG132 stock solution to 90 $\mu$ L 1 $\times$  proteasome assay buffer to give 100 $\mu$ L working solution

**Storage:** Stock and working solutions can be stored at -80°C for up to one month.

## AMC standard (optional)

AMC standard supplied as a solid (1.75 $\mu$ g)

- Prepare a 160 $\mu$ M AMC standard stock solution by dissolving solid in 62.5 $\mu$ L DMSO
- Prepare a 16 $\mu$ M AMC standard working solution by 10 $\times$  dilution in 1 $\times$  proteasome assay buffer – add 10 $\mu$ L 160 $\mu$ M AMC standard stock solution to 90 $\mu$ L 1 $\times$  proteasome assay buffer to give 100 $\mu$ L standard working solution

**Storage:** AMC standard solutions can be stored at -80°C for up to three months protected from light. Avoid multiple freeze-thaws

## End point assays (linear response)

In order to run meaningful end point assays, for example to screen for proteasome activators and inhibitors, the assay conditions used must give a linear response (RFU vs. time) over the course of the assay. When using the StressXpress®20S proteasome assay kit for the first time it is recommended that linear response for the starting point assay set-up given is verified for your system using the continuous analysis approach prior to application in end point assays.

## Proteasome activation by SDS

It has often been reported that SDS can be used to activate the 20S core proteasome. This should be treated with caution as results, upon attempted 20S proteasome activation by SDS, are at best mixed and can sometimes lead to inhibition of activity.

## Activators / inhibitors

Pre-incubate at desired concentration with proteasome in assay mixture at assay temperature for 15-30 minutes prior to initiation by addition of Suc-LLVY-AMC.

## AMC standard curve (optional)

Plotting an AMC standard curve allows proteasome activity assay response to be converted from relative fluorescence units (RFU) to defined units (e.g. pmol AMC) if required for subsequent kinetic data analysis.

## Control reactions

*Recommended:* Substrate only control: No 20S proteasome – adjust readings for background intrinsic substrate fluorescence and any auto-hydrolysis of Suc-LLVY-AMC. Positive control: 20S proteasome – demonstrate assay working / components active.

*Optional:* Inhibitor control: Test proteasome inhibition with MG132 control inhibitor – comparison for inhibitor assays. Vehicle control: Solvent for activator / inhibitor compounds (e.g. DMSO) – determine effect, if any, of carrier for test compounds on assay performance.

## Assay optimisation

Optimal assay conditions for specific experiments must be determined by the user. Adjustment of the following parameters may facilitate this process:

- 20S proteasome concentration, 2.5-20nM
- Suc-LLVY-AMC concentration, 50-500 $\mu$ M
- Assay incubation (end point) or monitoring (continuous) time, 15-60 minutes
- Assay temperature, increase from room temperature to 30-37°C
- Activator / inhibitor concentration and pre-incubation time

## ASSAY PROTOCOL

Key points for assay protocol:

- Protocol described for 96-well plate format, assays can be run and analysed in micro-cuvette format if required (up to 100x 50 $\mu$ L assays)
  - Recommended that all samples, controls and standards are run in triplicate
  - Ensure Suc-LLVY-AMC solution equilibrated to assay temperature before addition
1. Set-up the microplate reader for recording at 360nm excitation and 460nm emission.
  2. Prepare the required number of 20S proteasome assays (50 $\mu$ L final volume) in an appropriate number of wells by addition of assay components in the order listed in the table below
  3. Use the plate configuration table provided to aid in proper sample, control and standard identification
  4. Mix components gently and incubate plate at room temperature for 15-30 minutes
  5. Add Suc-LLVY-AMC substrate to initiate assays and mix reagents by shaking plate gently for 15 seconds


Component	20S proteasome	Substrate only (blank)	Control inhibitor
1x proteasome assay buffer*	40.0	45.0	35.0
25nM 20S proteasome	5.0	-	5.0
50 $\mu$ M MG132 (control inhibitor)	-	-	5.0
	Incubate for 15-30 minutes at room temperature		
1mM Suc-LLVY-AMC	5.0	5.0	5.0

\*Adjust proteasome assay buffer volume to allow for 20S proteasome, Suc-LLVY-AMC and activator / inhibitor volume used

6. Measure the fluorescence signal:
  - a. **For continuous kinetic reading:** Commence timecourse assay data acquisition IMMEDIATELY upon addition of Suc-LLVY-AMC substrate using microplate reader. Monitor RFU response for 15-60 minutes taking regular readings (see Figure 1 for control 20S proteasome assay)
  - b. **For end point reading:** Incubate the reaction for 15-60 minutes in the dark. Mix the reagents and measure fluorescence intensity using microplate reader

## Data Analysis

1. Following assay completion average the triplicate readings for each sample, control and substrate only (blank) assay and subtract the average blank value from each 20S proteasome assay data point to give RFU values
2. **For continuous kinetic analysis:**
  - a. Plot data as RFU vs. time for each sample / control assay. To convert RFU to concentration of AMC released first produce an AMC standard curve, as described in the following section
  - b. If non-linear determine the range of initial time points during which the assay response is linear
  - c. Obtain the slope of a line fit to the data, using an appropriate linear regression program, to give the initial reaction velocity ( $V_0$ ) in RFU / min
  - d. Perform appropriate data analyses as required, for example Michaelis-Menten kinetics to determine  $K_M$ ,  $k_{cat}$  and  $k_{cat}/K_M$  data or to calculate specific activity
3. **For end point analysis:**
  - a. Plot data as RFU versus concentration, for example to assess performance of test compounds.
  - b. Perform appropriate data analyses as required, for example determining % inhibition % or  $IC_{50}$



**Figure 1:** Time course of control 20S proteasome activity with Suc-LLVY-AMC substrate using the StressXpress® 20S proteasome assay kit – 20S proteasome (~2.5nM) was incubated with 100µM Suc-LLVY-AMC in proteasome assay buffer for 30 minutes at room temperature alongside a substrate only (blank) control reaction. Fluorescence measurements (RFU) were taken at 30 second intervals and plotted vs. time

## AMC standard curve (optional)


Perform an 'in-well' serial dilution of 16 $\mu$ M AMC standard working solution in proteasome assay buffer in triplicate as follows:

- Add 50 $\mu$ L proteasome assay buffer to wells A to H
- Add 50 $\mu$ L 16 $\mu$ M AMC standard working solution to well A
- Pipette 50 $\mu$ L from well A and mix it thoroughly with the contents of well B
- Pipette 50 $\mu$ L from well B and mix it thoroughly with the contents of well C etc. up to well G
- Discard the last 50 $\mu$ L pipetted from well G
- Leave well H as a blank

	AMC Standard ( $\mu$ M)		
	1	2	3
A	8.0	8.0	8.0
B	4.0	4.0	4.0
C	2.0	2.0	2.0
D	1.0	1.0	1.0
E	0.5	0.5	0.5
F	0.25	0.25	0.25
G	0.125	0.125	0.125
H	0.0	0.0	0.0

1. Measure fluorescence reading for AMC standards using microplate reader (360nm excitation, 460nm emission)
2. Average the triplicate readings for each AMC concentration
3. Plot RFU vs. AMC concentration to visualise results
4. Determine slope of the plot using an appropriate linear regression program





**Figure 2:** AMC Standard Curve – Representative plot of fluorescence measurements (RFU) of serial dilutions of AMC Standard at recommended concentrations

# RESOURCES

## Plate Configuration Table

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

## References

1. Saeki, Y. & Tanaka, K. *Methods in molecular biology* (Clifton, N.J.) 832, 315–37 (2012). PMID: 22350895
2. Frankland-Searby, S. & Bhaumik, S. R. *Biochimica et biophysica acta* 1825, 64–76 (2012). PMID: 22037302
3. Basler, M., Kirk, C. J. & Groettrup, M. *Current opinion in immunology* 1–7 (2012). PMID: 23219269
4. Löw, P. *General and comparative endocrinology* 172, 39–43 (2011). PMID: 21324320
5. Dennissen, F. J. a, Kholod, N. & Van Leeuwen, F. W. *Progress in neurobiology* 96, 190–207 (2012). PMID: 22270043
6. Geng, F., Wenzel, S. & Tansey, W. P. *Annual review of biochemistry* 81, 177–201 (2012). PMID: 22404630

## Warranty and Limitation of Remedy

StressMarq Biosciences Inc. makes **no warranty or guarantee** of any kind, whether written or oral, expressed or implied, including without limitation, any warranty of fitness for a particular purpose, suitability and merchantability, which extends beyond the description of the chemicals hereof. StressMarq **warrants only** to the original customer that the material will meet our specifications at the time of delivery. StressMarq will carry out its delivery obligations with due care and skill. Thus, in no event will StressMarq have **any obligation or liability**, whether in tort (including negligence) or in contract, for any direct, indirect, incidental or consequential damages, even if StressMarq is informed about their possible existence. This limitation of liability does not apply in the case of intentional acts or negligence of StressMarq, its directors or its employees.

Buyer's **exclusive remedy** and StressMarq's sole liability hereunder shall be limited to a refund of the purchase price, or at StressMarq's option, the replacement, at no cost to Buyer, of all material that does not meet our specifications.

Said refund or replacement is conditioned on Buyer giving written notice to StressMarq within thirty (30) days after arrival of the material at its destination. Failure of Buyer to give said notice within thirty (30) days shall constitute a waiver by Buyer of all claims hereunder with respect to said material.

**For further details, please refer to our Warranty and Refund Policy located on our website and in our catalog.**

## NOTES

This document is copyrighted. All rights are reserved. This document may not, in whole or part, be copied, photocopied, reproduced, translated, or reduced to any electronic medium or machine-readable form without prior consent, in writing, from StressMarq Biosciences Inc. ©08/28/2010, StressMarq Biosciences Inc., Victoria, BC Canada, All rights reserved. Printed in Canada.