

Anti-LAMP1

Catalog# SMC-140C/D

Size: 25/100µg

PO Box 30244, Suite 405,
3989 Quadra Street,
Victoria, BC V8X 5E1, Canada

This product is for *in vitro* research use only and is not intended for use in humans or animals

StressMarq

Biosciences Inc.

Orders • sales@stressmarq.com
Tel: • +1 250 294 9065
Fax: • +1 250 294 9025
Email • info@stressmarq.com
Web • www.stressmarq.com

Product	Mouse Anti-Lysosome-Associated Membrane Protein 1
Clone	Ly1C6
Immunogen	Rat liver lysosomal membrane preparations
Host and Subclass	Mouse monoclonal, IgG ₁
Cited Applications	WB (6), IP (5), ICC (1)
Specificity	Detects an ~120kDa protein corresponding to the molecular mass of LAMP1 on SDS-PAGE immunoblots
Species cross-reactivity	Rat, Hamster
Format	Mouse immunoglobulin in PBS pH 7.4, 0.09% azide, in 50% glycerol
Concentration and working dilution	1.0mg/mL, 1µg/ml was sufficient for detection by immunoblot analysis
Storage and stability	-20°C; 1 year+; shipped on cold packs or ambient

Scientific Background

Lysosome associated membrane proteins, or LAMP1 and LAMP2, are major constituents of the lysosomal membrane. The two have closely related structures, with 37% sequence homology (2). They are both transmembrane glycoproteins that are localized primarily in lysosomes and late endosomes. Newly synthesized molecules are mostly transported from the trans-Golgi network directly to endosomes and then to lysosomes. A second pathway involves the lamps being delivered from the Golgi to the cell surface, and then along the endocytic pathway to the lysosomes. A minor pathway involves transport via the plasma membrane (3).

Upon stimulation, a rapid translocation of intracellular LAMPs to the cell membrane is dependent on a carboxyl-terminal tyrosine based motif (YXXI) (1). If there is a disturbance in this spacing, lysosome localization of LAMP1 is abolished and the mutant protein then cycles between the membrane and the endosome (3). This stimulation has also been shown to have an associated

release of histamine, leukotriene C (4) and prostaglandin D (2), which shows that LAMP-1 and LAMP-2 are activation markers for normal mast cells (1). They have also been linked to the inflammatory response in that they promote adhesion of human peripheral blood mononuclear cells (PBMC) to vascular endothelium, and therefore possibly the adhesion of PBMC to the site of inflammation (4).

Selected References

1. Grutzkau, A. et al (2004) *Cytometry A*. 61(10): 62-68.
2. Furuta, K. et al. (1999) *EMBO J*. 17(5):1304-14.
3. Rohrer, J. et al. (1996) *J Cell Biol*. 132(4): 565-76.
4. Kannan, K., et al. (1996) *Cell Immunol*. 171: 10-19.
5. Lewis, V., et al. (1985) *J. Cell Biol*. 100: 1839-1847.
6. Jones, K.A., et al. (2004) *Exp Cell Res* 295(2): 512-524.

Certificate of Analysis

1.0 µg/ml was sufficient for detection of SMC-140 in rat liver microsomes by ECL immunoblot analysis

Material Safety Data Sheet

Anti-LAMP1 (Monoclonal Antibody) SMC-140

This product is for *in vitro* research use only and is not intended for use in humans or animals

The below information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. StressMarq shall not be held liable for any damage resulting from handling or from contact with the above product. See the Technical Specification, Packing Slip, Invoice, and Product Catalogue for additional terms and conditions of sale.

Hazardous Ingredients

The physical, chemical and toxicological properties of these components have not been fully investigated. It is recommended that all laboratory personnel follow standard laboratory safety procedures when handling this product. Safety procedures should include wearing OSHA approved safety glasses, gloves and protective clothing. Direct physical contact with this product should be avoided.

<u>Known Hazardous Components</u>	<u>CAS Number</u>	<u>Percent</u>
Sodium Azide	26628-22-8	0.09

Physical Data

This product consists of mouse immunoglobulin in TBS containing 0.09% azide in 50% glycerol shipped on gel packs. The physical properties of this product have not been investigated thoroughly.

Fire and Explosion Hazard and Reactivity Data

NOT APPLICABLE

Toxicological Properties

May be harmful by inhalation, ingestion, or skin absorption. The toxicological properties of this product have not been investigated thoroughly. Exercise due caution.

Preventative Measures

Wear chemical safety goggles and compatible chemical-resistant gloves. Avoid inhalation, contact with eyes, skin or clothing.

Spill and Leak Procedures

Observe all federal, state and local environmental regulations.

- Wear protective equipment.
- Absorb on sand or vermiculite and place in closed containers for disposal.
- Dispose or mix the material with a combustible solvent and burn in a chemical incinerator equipped with an afterburner and scrubber.

First Aid Measures

- If swallowed, wash out mouth with water, provided person is conscious. Call a physician.
- In case of skin contact, flush with copious amounts of water for at least 15 minutes. Remove contaminated clothing and shoes. If a rash or other irritation develops, call a physician.
- If inhaled, remove to fresh air. If breathing becomes difficult, call a physician.
- In case of eye contact, flush with copious amounts of water for at least 15 minutes while separating the eyelids with fingers. Call a physician.

Authorized: StressMarq Biosciences Inc.
Creation Date: 11/30/07