

Anti-GRP78 (BiP)

Catalog# SPC-107A/B
Size: 50/200µl

PO Box 30244, Suite 405,
3989 Quadra Street,
Victoria, BC V8X 5E1, Canada

This product is for *in vitro* research use only and is not intended for use in humans or animals

StressMarq

Biosciences Inc.

Orders • sales@stressmarq.com
Tel: • +1 250 294 9065
Fax: • +1 250 294 9025
Email • info@stressmarq.com
Web • www.stressmarq.com

Product	Rabbit anti-GRP78 (BiP) antibody; polyclonal
Clone	N/A
Immunogen	Rat GRP78 (BiP) synthetic peptide conjugated to KLH
Host and Subclass	Rabbit
Cited Applications	WB
Specificity	Detects a protein with a mass of ~78 kDa corresponding to Grp78
Species Cross-Reactivity	Human, Mouse, Rat, Cow, Fungus (<i>Cryptosporidia parvifolia</i>), Hamster, Monkey, Rabbit, <i>Xenopus</i>
Format	Rabbit antiserum
Working dilution	Dilution of 1:2000 for WB using ECL.
Storage and stability	-20°C; 1 year+; shipped on cold packs

Scientific Background

GRP78 is a ubiquitously expressed, 78-kDa glucose-regulated protein, and is commonly referred to as an immunoglobulin chain binding protein (BiP). The BiP proteins are categorized as stress response proteins because they play an important role in the proper folding and assembly of nascent protein and in the scavenging of misfolded proteins in the endoplasmic reticulum lumen. Translation of BiP is directed by an internal ribosomal entry site (IRES) in the 5' nontranslated region of the BiP mRNA. BiP IRES activity increases when cells are heat stressed (1).

GRP78 is also critical for maintenance of cell homeostasis and the prevention of apoptosis (2). Luo et al. have provided findings that suggest GRP78 is essential for embryonic cell growth and pluripotent cell survival (3). In terms of diseases, GRP78 has been shown to be a reliable biomarker of hypoglycemia, to serve a

neuroprotective function in neurons exposed to glutamate and oxidative stress (4), and its protein levels are reduced in the brains of Alzheimer's patients (5). Also, the induction of the GRP78 protein that results in severe glucose and oxygen deprivation could possible lead to drug resistance to anti-tumor drugs (6, 7).

Selected References

1. Cho, S. *et al.* (2007). *Mol Cell Biol* 27(1): 368-83.
2. Yang, Y. *et al.* (1998) *J Biol Chem* 273: 25552-25555.
3. Luo, S. *et al.* (2006) 26 (15): 5688-97.
4. Yu, Z. *et al.* (1999) *Exp Neurol*. 15: 302-314.
5. Koomagi, R. *et al.* (1999) *Anticancer Res*. 19: 4333-4336.
6. Laquerre, S. *et al.* (1998) *J. Virology* 72: 4940-4949.
7. Dong, D. *et al.* (2005) *Cancer Res* 65(13): 5785-91.

Certificate of Analysis

0.5 µg/ml of SPC-107 was sufficient for detection of Grp78 in 10ug of rat tissue lysate by ECL immunoblot analysis.

Material Safety Data Sheet

GRP78 (BIP) Polyclonal SPC-107

This product is for *in vitro* research use only and is not intended for use in humans or animals

The below information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. StressMarq shall not be held liable for any damage resulting from handling or from contact with the above product. See the Technical Specification, Packing Slip, Invoice, and Product Catalogue for additional terms and conditions of sale.

Hazardous Ingredients

The physical, chemical and toxicological properties of these components have not been fully investigated. It is recommended that all laboratory personnel follow standard laboratory safety procedures when handling this product. Safety procedures should include wearing OSHA approved safety glasses, gloves and protective clothing. Direct physical contact with this product should be avoided.

<u>Known Hazardous Components</u>	<u>CAS Number</u>	<u>Percent</u>
-----------------------------------	-------------------	----------------

Physical Data

This product consists of rabbit antiserum shipped on gel packs. The physical properties of this product have not been investigated thoroughly.

Fire and Explosion Hazard and Reactivity Data

NOT APPLICABLE

Toxicological Properties

May be harmful by inhalation, ingestion, or skin absorption. The toxicological properties of this product have not been investigated thoroughly. Exercise due caution.

Preventative Measures

Wear chemical safety goggles and compatible chemical-resistant gloves. Avoid inhalation, contact with eyes, skin or clothing.

Spill and Leak Procedures

Observe all federal, state and local environmental regulations.

- Wear protective equipment.
- Absorb on sand or vermiculite and place in closed containers for disposal.
- Dispose or mix the material with a combustible solvent and burn in a chemical incinerator equipped with an afterburner and scrubber.

First Aid Measures

- If swallowed, wash out mouth with water, provided person is conscious. Call a physician.
- In case of skin contact, flush with copious amounts of water for at least 15 minutes. Remove contaminated clothing and shoes. If a rash or other irritation develops, call a physician.
- If inhaled, remove to fresh air. If breathing becomes difficult, call a physician.
- In case of eye contact, flush with copious amounts of water for at least 15 minutes while separating the eyelids with fingers. Call a physician.

Authorized: StressMarq Biosciences Inc.
Creation Date: 12/12/2008